


Marquette-Hellige CardioSys Cycle Interface

You can setup the Cycle with your stress system using the following directions. You will need the RS-232 cable specified below. It establishes communication between the monitor and your stress system. This cable allows your stress system to prompt the monitor when it needs a BP measurement. In addition, it allows the BP measurements taken by the Cycle to be transferred to your stress system display and reports.



RS-232 cable
Part# 91-0013-00



1. Plugging in the RS-232 cable

Connect: the RS-232 cable, SunTech part #91-0013-00 (9 pin female to 9 pin female)
From: the RS-232 connection on the connector panel of the Cycle (9 pin female)
To: the COM 1, 2, or 3 port on the back of the stress system (9 pin female)

2. Setting up the Cycle monitor

- From the measurement view, press and hold the **SELECT** and **ADVANCE** buttons simultaneously for 3 seconds. The System Menu will be displayed.
- Press the **MENU** button once to proceed to the Communication Screen. Press **ADVANCE** until *Cardiosys* is displayed.
- Press **MENU** twice to return to the measurement view.

3. Setting up the stress system

- Turn the CardioSys on and at the main screen, click on the “running man” icon. This will present the **Retrieve Patient** screen.
- Select any patient and press **OK**. At the **Test Demographics** screen press **Cancel**.
- Select **Defaults** and proceed to **Page 2**. Select the **1mv signal** by ensuring that an ‘x’ appears in the box.
- Select the **BP Monitor** as **BP METER 1**
- For **Port**, choose the COM port you have connected the cable to i.e. **COM 1**, **COM 2** or **COM 3**
- Select **OK** to exit.