ENGLISH

SunTech Medical's Clinical Grade Orbit-K Blood Pressure Cuffs

Index Line.

Cuff index line must fall within range markings.

Artery symbol and arrow should be placed over brachial or femoral artery.

Symbol indicating product contains no PVC (polyvinyl chloride) material.

Symbol indicating product contains no latex material.

Product in compliance with Council Directive 93/42/EEC Medical Devices Directive.

Symbol indicating arm circumference.

Symbol indicating to refer to "Instructions for Use"

RANGES

Size:	PN:	Range:
Sm Adult	98-0061-01	18-27 cm
Adult	98-0061-02	25-35 cm
Adult Plus	98-0061-05	27-40 cm
Lg Adult	98-0061-03	32-44 cm

It is important that the cuff is properly fitted to the patient's arm, and that the microphone is placed over the brachial artery (on the inside of the arm between the bicep and tricep). Improper cuff sizing and a misplaced microphone can lead to missed or poor readings and accuracy. Please refer to the Tango+ and Tango M2 User Manual for further information about this cuff that is used exclusively with the Tango Family Product Line.

POSSIBLE ADVERSE REACTIONS

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechiae) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

WARNINGS & CONTRAINDICATIONS

Federal (US) law restricts this device to sale by or on the order of a physician.

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through one of the hose openings in the cuff.

Cuffs are intended for use with non-invasive blood pressure monitors only.

DO NOT apply the BP cuff to a limb being used for IV infusions as the cuff inflation can block the infusion, causing harm to the patient.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

WARNING: Regarding the effect of blood flow interference and resulting harmful injury to the PATIENT caused by continuous CUFF pressure due to connection tubing kinking.

DO NOT apply the BP cuff over a wound as this can cause further injury.

DO NOT apply the BP cuff to the arm on the side of a single mastectomy. In the case of double mastectomy use the side of the least dominate arm.

APPLICATION

Check that the cuff is the correct size:

1. Fold the grey sleeve inside the blue cuff (away from the Velcro strip).
2. Wrap the cuff around the patient's upper arm.

3. Make sure the INDEX (the end of the cuff) falls

within the RANGE (printed inside the cuff).

4. If the INDEX falls outside the RANGE, select a new cuff size.

WARNING: Using an incorrect cuff size could result in erroneous and misleading BP measurements!

PROPER CUFF & MICROPHONE PLACEMENT

1. Locate the brachial artery, between the bicep and the tricep of the upper arm. The left arm is preferred.

2. Slide the cuff sleeve up the patient's arm, with the "ARTERY" marker pointing down the arm.

3. There is a microphone located under the "ARTERY" marker. Make sure the microphone is placed on the inner portion of the arm, directly over the brachial artery between the edge of the cuff and the elbow.

4. Wrap the cuff around the arm and secure. Use the wrist straps to secure the cables to the patient's wrist.

NOTE: It is important that the cuff is properly fitted to the patient's arm, and that the microphone is placed over the brachial artery in between the bicep and tricep on the inner portion of the upper arm.

CLEANING INSTRUCTIONS

NOTE: Orbit-K cuff should be cleaned at the conclusion of each stress test.

1. Use a medical grade mild disinfectant wipe on the cuff or spray a cleaning solution onto a cloth and wipe the cuff. Afterwards, lay flat or line dry.

2. Remove the bladder and microphone from the outer shell of the Orbit-K cuff. Machine wash the shell in warm water with a mild detergent (50-140°F or 10-60°C). Lay flat or line dry the cuff. Do not place the cuff in a dryer.

CAUTION: Do not machine wash bladder or microphone.

The bladder needs to be inserted back into the cuff sleeve so the pneumatic hose portion of the bladder is outside the sleeve. Please note that the pneumatic hose connection should face downward when using the Orbit-K cuff on either the right or the left arm.

WARRANTY

SunTech Medical, Inc. provides to the original purchaser the following limited warranty from date of invoice.

Orbit-K Cuffs: 6 months

Accessories, i.e. patient cables, microphone, disposables: 90 days

This limited warranty is void if the instrument has been damaged by accident, misuse, negligence, act of God or serviced by any person not authorized by SunTech Medical, Inc.

It is recommended that you replace the Orbit-K cuffs, microphones and patient cable annually to maintain measurement accuracy.

Cuffs are intended for use with non-invasive blood pressure monitors only.

DO NOT apply the BP cuff to a limb being used for IV infusions as the cuff inflation can block the infusion, causing harm to the patient.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

WARNING: Regarding the effect of blood flow interference and resulting harmful injury to the PATIENT caused by continuous CUFF pressure due to connection tubing kinking.

DO NOT apply the BP cuff over a wound as this can cause further injury.

DO NOT apply the BP cuff to the arm on the side of a single mastectomy. In the case of double mastectomy use the side of the least dominate arm.

APPLICATION

Check that the cuff is the correct size:

1. Fold the grey sleeve inside the blue cuff (away from the Velcro strip).

2. Wrap the cuff around the patient's upper arm.

3. Make sure the INDEX (the end of the cuff) falls

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SunTech Medical, Inc. offre à l'acheteur d'origine la garantie limitée suivante à compter de la date de la facture.
SunTech Medical Ltd. Europe