All Purpose Durable BP Cuff

English

English

Suntech Medical's All Purpose Clinical Grade Blood Pressure Cuffs

	Index Line
	Cuff index line must fall within range markings
	Artery symbol and arrow should be placed over brachial or femoral artery
	Not made with PVC
	Not made with natural rubber latex
	Product in compliance with Council Directive 93/42/EEC Medical Devices Directive
	Symbol indicating arm circumference
	Symbol indicating manufacturer
	Symbol indicating lot code of cuff
	Symbol indicating Caution

RANGES / COLORS

Size:	Color:	Range:
Infant	Orange	8-13 cm
Child	Green	12-19 cm
Child Long	Green	12-19 cm
Sm Adult	Royal Blue	17-25 cm
Sm Adult Long	Royal Blue	17-25 cm
Adult	Navy Blue	23-33 cm
Adult Long	Navy Blue	23-33 cm
Lg Adult	Burgundy	31-40 cm
Lg Adult Long	Burgundy	31-40 cm
Thigh	Brown	38-50 cm

WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

The cuff should not be placed on the patient's arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

CAUTIONS

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through one of the hose openings in the cuff.

Minimize limb movement during the measurement.

The cuff should not be applied over a wound as this can cause further injury

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy

Do not machine wash the cuff bladder. Water could be trapped in the cuff and cause damage to the NIBP module and/or inaccurate BP readings Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress. A compressed or kinked connection hose may

cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient

Using an incorrect cuff size could result in erroneous and misleading BP measurement results

APPLICATION

Follow the application instructions for use to ensure the correct size cuff for the patient. Failure to do so will adversely affect the accuracy of the reading.

- Place open cuff around the inner portion of the upper arm (or thigh).
- Align artery symbol ARTERY to the brachial (or femoral) artery.
- Use the RANGE indicator with the INDEX line to check that the arm falls within the specified range of that cuff. If it does not, select a cuff that better accommodates the limb circumference.
- Wrap the cuff snugly around the arm (or thigh).

Promptly remove cuff from patient when monitoring is not in progress.

CLEANING INSTRUCTIONS

The following cleaning methods have been applied 20 times to the cuff without any apparent negative effects.

The cuff may be sprayed with a mild disinfectant solution (e.g. Cidezyme® ENZOL®, or 10% bleach solution), rinsed with distilled water and line dry. Ensure that no liquid enters tubing.

OR

To machine wash the cuff, remove the bladder and fully engage the hook and the loop. Machine wash warm with a mild detergent (50 – 130° F or 1 - 54° C) and line dry.

Disposal: Please dispose of cuff according to local regulations

ENVIRONMENTAL CONDITIONS

Operating Ranges
Temp: 0 - 50°C
RH: 15 – 95% non-condensing
Storage Ranges
Temp: -20 - 65°C
RH: 15 – 95% non-condensing

WARRANTY

SunTech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. This limited warranty covers the no charge replacement of the cuff under normal wear and tear conditions when returned to the attention of Service Department at the address below depending on location. Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Service Department at of the addresses below.

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Suntech Medical, Ltd. Europe