One Piece

Durable BP Cuff

English

Suntech Medical's One Piece Clinical Grade Blood Pressure Cuffs

-  Index Line
-  Cuff index line must fall within range markings
-  Artery symbol and arrow should be placed over brachial or femoral artery
-  Not made with PVC
-  Not made with natural rubber latex
-  Product in compliance with Council Directive 93/42/EEC Medical Devices Directive
-  Symbol indicating arm circumference
-  Symbol indicating lot code of cuff
-  Symbol indicating manufacturer
-  Symbol indicating Caution

RANGES / COLORS

Size:	Color:	Range:
Infant	Orange	8-13 cm
Child	Green	12-19 cm
Sm Adult	Royal Blue	17-25 cm
Adult	Navy Blue	23-33 cm
Adult Long	Navy Blue	23-33 cm
Adult Plus	Navy Blue	28-40 cm
Lg Adult	Burgundy	31-40 cm
Lg Adult Long	Burgundy	31-40 cm
Lg Adult Plus	Burgundy	40-55 cm
Thigh	Brown	38-50 cm

WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb. The cuff should not be placed on the patient's arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm. If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

CAUTIONS

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner. The cuff should not be applied over a wound as this can cause further injury.

Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy.

Minimize limb movement during the measurement. Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress.

A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.

Using an incorrect cuff size could result in erroneous and misleading BP measurement results.

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SunTech Medical, Ltd.
SunTech Medical, Inc.