ENGLISH

Suntech Medical's All Purpose Clinical Grade Blood Pressure Cuffs

- Index Line.
- Cuff index line must fall within range markings.
- Artery symbol and arrow should be placed over brachial or femoral artery.
- Symbol indicating product contains no PVC (polyvinyl chloride) material.
- Symbol indicating product contains no latex material.
- Product in compliance with Council Directive 93/42/EEC Medical Devices Directive.
- Symbol indicating arm circumference

RANGES / COLORES

Size:	Color:	Range:
Thigh	Brown	38-50 cm
Lg Adult Long	Burgundy	31-40 cm
Lg Adult	Burgundy	31-40 cm
Adult Long	Navy Blue	23-33 cm
Adult	Navy Blue	23-33 cm
Sm Adult Long	Royal Blue	17-25 cm
Sm Adult	Royal Blue	17-25 cm
Child Long	Green	12-19 cm
Child	Green	12-19 cm
Infant	Orange	8-13 cm

WARNINGS & INDICATIONS
 Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through one of the hose openings in the cuff.

Replace cuff when worn or damaged.

Do not machine wash bladder.

If patient has a mastectomy, place the cuff on the arm opposite the side of the mastectomy. If a patient has a double mastectomy, the cuff must be placed on the least dominate arm.

If a patient is undergoing intravascular access or therapy, or an arterio-venous (A-V) shunt, due to the presence of temporary interference of blood flow, do not apply the cuff on that arm. Potential injury to patient may occur.

The use of cuff while using other monitoring medical equipment could affect the accuracy or operation of said equipment.

Federal (US) law restricts this device to sale by or on the order of a physician.

Cuffs are intended for use with non-invasive blood pressure monitors only.

Follow the application instructions for use to ensure the correct size cuff for the patient. Failure to do so will adversely affect the accuracy of the reading.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Minimize limb movement during the measurement.

Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress.

- Do not compress or pinch cuff tubing.
- Ensure patient's skin is intact before applying the cuff.
- Promptly remove cuff from patient when monitoring is not in progress.
- During frequent or extended measurement, check cuff site and limb to ensure proper blood flow. Do not apply cuff where circulation may be compromised.
- APPLICATION**
 - Wrap the cuff around the patient's upper arm midway between the elbow and the shoulder.
 - Ensure the ARTERY arrow is over the brachial artery, between the biceps and triceps muscles on the inside of the arm.
 - Use the RANGE indicator with the INDEX line to check that the arm falls within the specified range for that cuff. Use appropriate larger or smaller cuff if necessary.
 - Wrap the cuff snugly around the arm.

CLEANING INSTRUCTIONS
 The following cleaning methods have been applied 20 times to the cuff without any apparent negative effects.

The cuff may be sprayed with a mild disinfectant solution (e.g. Cidezyme® ENZOL®, or 10% bleach solution), rinsed with distilled water and line dry. Ensure that no liquid enters tubing.

OR

To machine wash the cuff, remove the bladder and fully engage the hook and the loop. Machine wash warm with a mild detergent (50 – 130° F or 1 - 54° C) and line dry.

ENVIRONMENTAL CONDITIONS
 Operating Ranges

- Temp: 0 - 50°C
- RH: 15 – 95% non-condensing
- Storage Ranges
- Temp: -20 - 65°C
- RH: 15 – 95% non-condensing

WARRANTY
 SunTech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. This limited warranty covers the no charge replacement of the cuff under normal wear and tear conditions when returned to the attention of Service Department at the address below depending on location.

Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Service Department at the addresses below.

www.SunTechMed.com
 SunTech Medical, Inc.
 507 Airport Boulevard, Suite 117
 Morrisville, NC 27560-8200, USA
 TEL: 1.800.421.8626 or 1.919.654.2300
 FAX: 1.919.654.2301

SunTech Medical, Ltd. Europe
 Oakfield Industrial Estate
 Stanton Harcourt Road
 Eynsham,
 Oxfordshire OX29 4TS
 England
 TEL: 44.1 865.884 234
 FAX: 44.1 865.884 235

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SunTech Medical, (Shenzhen) Co., Ltd.