Soft Disposable

Single Patient Use BP Cuff

English

Suntech Medical's Disposable Clinical Grade Blood Pressure Cuffs

Index Line

Cuff index line must fall within range markings

Artery symbol and arrow should be placed over brachial or femoral artery

Not made with natural rubber latex

Product in compliance with Council Directive 93/42/EEC Medical Devices Directive

Symbol indicating arm circumference

Symbol indicating lot code of cuff

Symbol indicating manufacturer

Symbol indicating Caution

RANGES / COLORS

| Size: | Color: | Range: |
|---------------|-----------|----------|
| Infant | Orange | 8-13 cm |
| Child | Green | 12-19 cm |
| Child Long | Green | 12-19 cm |
| Sm Adult | Turquoise | 17-25 cm |
| Sm Adult Long | Turquoise | 17-25 cm |
| Adult | Navy Blue | 23-33 cm |
| Adult Long | Navy Blue | 23-33 cm |
| Adult Plus | Navy Blue | 28-40 cm |
| Lg Adult | Burgundy | 31-40 cm |
| Lg Adult Long | Burgundy | 31-40 cm |
| Lg Adult Plus | Burgundy | 40-55 cm |
| Thigh | Brown | 38-50 cm |

Neonatal:

| Size: | Color: | Range: |
|-------|-----------|---------|
| 1 | Orange | 3-6 cm |
| 2 | Turquoise | 4-8 cm |
| 3 | Green | 6-11 cm |
| 4 | Navy Blue | 7-13 cm |
| 5 | Burgundy | 8-15 cm |

WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

The cuff should not be placed on the patient's arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

CAUTIONS

Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy

Minimize limb movement during the measurement.

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner

The cuff should not be applied over a wound as this can cause further injury

Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress.

A compressed or kinked connection hose may

cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient

Using an incorrect cuff size could result in erroneous and misleading BP measurement results

APPLICATION

Follow the application instructions for use to ensure the correct size cuff for the patient. Failure to do so will adversely affect the accuracy of the reading.

- Place open cuff around the inner portion of the upper arm (or thigh).
- Align artery symbol ARTERY to the brachial (or femoral) artery.
- Use the RANGE indicator with the INDEX line to check that the limb falls within the specified range of that cuff. If it does not, select a cuff that better accommodates the limb circumference.
- Wrap the cuff snugly around the arm (or thigh).

Promptly remove cuff from patient when monitoring is not in progress.

Disposal: Please dispose of cuff according to local regulations

ENVIRONMENTAL CONDITIONS

Operating Ranges
Temp: 0 - 50°C
RH: 15 – 95% non-condensing
Storage Ranges
Temp: -20 - 65°C
RH: 15 – 95% non-condensing

WARRANTY

SunTech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. The warranty shall only cover out-of-box failures, as this is a Single Use product. It is not intended for repeated use. Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Service Department at the addresses below.

SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

SunTech Medical, Ltd. Europe
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England

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