


Quinton Q5000 Tango Interface Notes

You can setup the Tango with your stress system using the following directions. You will need the ECG trigger and RS-232 cables specified below.

- The ECG trigger cable is required for the Tango monitor to be able to take a BP measurement.
- You will need the RS-232 cable to establish communication between the Tango and your stress system. This cable allows your stress system to prompt the Tango when it needs a BP measurement. In addition, it allows the BP measurements taken by the Tango to be transferred to your stress system display and reports.



1. Plugging in the ECG trigger cable
Connect: the ECG trigger cable, SunTech part #91-0018-00 (1/8" phono plug to BNC)
From: the BNC External ECG connection on the rear panel of the Tango (BNC)
To: the QRS SYNC connection on the back of the stress system (1/8" phono plug); if this connection is occupied, then use either
 - The connection labeled CHANNEL ONE OUT or
 - A dual mini phone jack adapter to allow two connections to the QRS SYNC connection.
2. Plugging in the RS-232 cable
Connect: the RS-232 cable, SunTech part #91-0020-00 (9 pin female to 9 pin male)
From: the RS-232 connection on the rear panel of the Tango (9 pin female)
To: the AUX 1 connection on the back of the stress system (9 pin male)

3. Setting up the Tango monitor
 - a. Turn on the Tango using the switch at the rear top left corner of the Tango.
 - b. Once the operating screen is displayed, press the *MENU/SELECT* button. Use the UP or DOWN buttons to select **Utilities** and press the menu button.
 - c. Use the arrow buttons to select **Device** and press the *MENU/SELECT* button to edit. Choose **QUINTON** using the arrow buttons and press the select button to confirm the choice.
 - d. Select **ECG trigger** and press the select button to edit. Choose **DIGITAL**↑ and press the select button to confirm the choice. If you are connected to **CHANNEL ONE OUT**, select **ANALOG** for the ECG trigger.
 - e. Select **Exit** twice to return to the operating screen.
4. Setting up the stress system
 - a. To setup the Tango as the BP monitor, choose **SET UP** at the master menu, then **SERIAL PORT CONFIGURATION**.
 - b. Select AUX 1/BP with **410/412 BP MONITOR** as the set up.
 - c. Select **PREVIOUS MENU**, then **MASTER MENU** to go to main screen.
 - d. To program the stress system when to take BP readings, choose **SELECT PROCEDURE** from the master menu.
 - e. Select **EDIT PROCEDURE** twice in the upper right corner and then the procedure you wish to edit.
 - f. Select **NEXT MENU** twice, scroll through the menu to **SELECT TIMED EXERCISE EVENTS**. Use the up and down arrows to select the time intervals at which you would like to take BP readings, or type in actual times.
 - g. Select **BP** for adding the event.
 - h. Select **NEXT MENU**, then **RECOVERY EVENTS**, and follow steps g and h to setup similarly.
 - i. Select **SAVE**, then go to **RECOVERY MENU**, and then to **MASTER MENU** to return.
 - j. Start a test or use the keyboard command to initiate a BP reading using the Q5000 to complete the setup.