Vinyl

Short-term Use BP Cuff

English

Warnings

Suntech Medical’s Vinyl Clinical Grade Blood Pressure Cuffs

	Index Line
	Cuff index line must fall within range markings

 Artery symbol and arrow should be placed over brachial or femoral artery

 Not made with natural rubber latex

 Product in compliance with Council Directive 93/42/EEC Medical Devices Directive

 Symbol indicating arm circumference

 Symbol indicating lot code of cuff

 Symbol indicating manufacturer

 Symbol indicating Caution

 Indicates device sale only by or on the order of a licensed practitioner

 Indicates keep dry

RANGES / COLORS		
Size:	Color:	Range:
Infant	Orange	8-13 cm
Child	Green	12-19 cm
Sm Adult	Royal Blue	17-25 cm
Sm Adult Plus	Royal Blue	18-29 cm
Adult	Navy Blue	23-33 cm
Adult Long	Navy Blue	23-33 cm
Adult Plus	Navy Blue	28-40 cm
Lg Adult	Burgundy	31-40 cm
Lg Adult Long	Burgundy	31-40 cm
Lg Adult Plus	Burgundy	40-55 cm
Thigh	Brown	38-50 cm

WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

The cuff should not be placed on the patient’s arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

CAUTIONS

Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy

Minimize limb movement during the measurement.

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner

The cuff should not be applied over a wound as this can cause further injury

Avoid contact with the cuff, other than that of the patient’s limb, while measurement is in progress.

A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood

Warnings

flow interference and potentially harmful injury to the patient

Using an incorrect cuff size could result in erroneous and misleading BP measurement results

APPLICATION
Follow the application instructions for use to ensure the correct size cuff for the patient. Failure to do so will adversely affect the accuracy of the reading.

- Place open cuff around the inner portion of the upper arm (or thigh).
- Align artery symbol **ARTERY** to the brachial (or femoral) artery.
- Use the **RANGE** indicator with the **INDEX** line to check that the limb falls within the specified range of that cuff. If it does not, select a cuff that better accommodates the limb circumference.
- Wrap the cuff snugly around the arm (or thigh).

Promptly remove cuff from patient when monitoring is not in progress.


CLEANING INSTRUCTIONS

The following cleaning methods have been applied 20 times to the cuff without any apparent negative effects.

The cuff may be sprayed with a mild disinfectant solution (e.g. Cidezyme®, ENZOL®, or 10% bleach solution), rinsed with distilled water and line dry. Ensure that no liquid enters tubing. OR The cuff may be wiped down with a mild disinfectant wipe (e.g. Sani-Cloth®) and line dry.

Disposal: Please dispose of cuff according to local regulations

ENVIRONMENTAL CONDITIONS

	Operating Ranges
	Temp: 0 - 50°C
	RH: 15 – 95% non-condensing
	Storage Ranges
	Temp: -20 - 65°C
	RH: 15 – 95% non-condensing

WARRANTY

SunTech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. The warranty shall only cover out-of-box failures, as this is a limited use product. It is not intended for extended use. Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Ser-vice Department at of the addresses below.

SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

SunTech Medical, Ltd. Europe
Oakfield Industrial Estate 
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235

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