

ENGLISH

Suntech Medical's Clinical Grade Ambulatory Blood Pressure Cuffs


Index Line.


←Range→ Cuff index line must fall within range markings.


Artery symbol and arrow should be placed over brachial or femoral artery.


Symbol indicating product contains no latex material.

Product in compliance with Council Directive 93/42/EEC Medical Devices Directive.

Symbol indicating arm circumference

Symbol indicating to refer to "Instructions for Use"

RANGES	PN:	Range:
1	98-0239-01	18-28 cm
2	98-0239-02	26-34 cm
3	98-0239-03	32-44 cm
4	98-0239-04	42-55 cm

It is important that the cuff is properly fitted to the patient's arm, and that the artery symbol on the cuff is placed over the brachial artery (on the inside of the arm between the bicep and tricep)! Improper cuff sizing or placement can lead to missed or poor readings and accuracy. Please refer to the Oscar 2 ABPM System User Manual for further information about this cuff that is used exclusively with the Oscar Family Product Line.

POSSIBLE ADVERSE REACTIONS
Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

WARNINGS & CONTRAINDICATIONS
Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

WARNING: Do not apply the BP cuff over a wound as this can cause further injury.

WARNING: Do not apply the BP cuff to the arm on the side of a single mastectomy. In the case of double mastectomy use the side of the least dominate arm.

WARNING: Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.

WARNING: Do not attach the BP cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

CAUTION: A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.

Do not use the BP cuff if there are any signs of damage. Failure to do so could affect measurement accuracy.

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through the hose opening in the cuff.

Cuffs are intended for use with non-invasive blood pressure monitors only.

APPLICATION
Check that the cuff is the correct size:

1. Fold the grey sleeve inside the blue cuff (away from the Velcro strip).
2. Wrap the cuff around the patient's upper arm.
3. Make sure the INDEX (the end of the cuff) falls within the RANGE (printed inside the cuff).
4. If the INDEX falls outside the RANGE, select a new cuff size.

CAUTION: Using an incorrect cuff size could result in erroneous and misleading BP measurements!

PROPER CUFF PLACEMENT

1. Locate the brachial artery, between the bicep and the tricep of the upper arm. The left arm is preferred.
2. Slide the cuff sleeve up the patient's arm, with the "ARTERY" marker pointing down the arm.
3. Make sure the "ARTERY" marker is placed on the inner portion of the arm, directly over the brachial artery between the edge of the cuff and the elbow.
4. Wrap the cuff around the arm and secure.

NOTE: It is important that the cuff is properly fitted to the patient's arm, and is placed over the brachial artery in between the bicep and tricep on the inner portion of the upper arm.

CLEANING INSTRUCTIONS
NOTE: SunTech Medical ABPM cuffs should be cleaned at the conclusion of each use or study.

Remove the bladder from the outer shell of the SunTech Medical ABPM cuff. Machine wash the shell in warm water with a mild detergent (50-140°F or 10-60°C). Lay flat or line dry the cuff. Do not place the cuff in a dryer.

CAUTION: Do not machine wash bladder of the cuff. The bladder needs to be inserted back into the cuff sleeve so the pneumatic hose portion of the bladder is outside the sleeve. Please note that the pneumatic hose connection should face upward when using the SunTech Medical ABPM cuff on either the right or the left arm.

Please dispose of cuff according to local regulations.

Operating Ranges	Storage Ranges
Temp: 0-50° C	Temp: -20 - 70° C
RH: 15 - 95%	RH: 15 - 95%

WARRANTY
SunTech Medical, Inc. provides to the original purchaser the following limited warranty from date of invoice.

SunTech Medical ABPM Cuffs: 6 months
Accessories, i.e. pouches, straps, disposable: 90 days
This limited warranty is void if the instrument has been damaged by accident, misuse, negligence, act of God or serviced by any person not authorized by SunTech Medical, Inc.

It is recommended that you replace the SunTech Medical ABPM cuffs annually to maintain measurement accuracy.

www.SunTechMed.com
SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200, USA
TEL: 1.800.421.8626 or 1.919.654.2300
FAX: 1.919.654.2301

SunTech Medical, Ltd. Europe
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
TEL: 44.1.865.884.234
FAX: 44.1.865.884.235

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SunTech Medical, (Shenzhen) Co., Ltd.