Eclipse

Durable BP Cuff
English

Symbol	Definitions	Standard/Source
	Artery symbol and arrow should be placed over brachial or femoral artery	Manufacturer
	Not made with PVC	Manufacturer
	Not made with natural rubber latex	Manufacturer
	Product in compliance with (EU)2017/745 MDR	EU Directive
	Symbol indicating arm circumference	Manufacturer
	Symbol indicating manufacturer	ISO 7000-3082
	Symbol indicating lot code of cuff	ISO 15223 3.13
	Symbol indicating Caution	ISO 7000-0434A
	Warning message	ISO 7010-W001
	Medical Device	Medical Device Regulations
	Device sale only by or on the order of a licensed practitioner	FDA
	Indicates the temperature limits to which the medical device can be safely exposed	ISO 7000-0632
	Indicates the range of humidity to which the medical device can be safely exposed	ISO 7000-2620
	To identify the country of manufacture of products. The "CC" shall be replaced by the two-letter country code defined in ISO3166-1.	IEC 60417-6049
	Indicates the model number or type number of a product	IEC 60417-6050

RANGES	
Pediatric	6-22 cm
Small Adult	21-29 cm
Adult	28-37 cm
Large Adult	36-46 cm

WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

The cuff should not be placed on the patient's arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm. If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury. Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient. It is recommended that the user monitor the patient for signs of prolonged impairment of circulation of the blood of the patient.

The cuff should not be applied over a wound as this can cause further injury.

Performing too frequent measurements can cause injury to the PATIENT due to blood flow interference.

CAUTIONS

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through one of the hose openings in the cuff. Minimize limb movement during the measurement.

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner

Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy

Do not machine wash the cuff bladder. Water could be trapped in the cuff and cause damage to the NIBP module and/or inaccurate BP readings

Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress.

Using an incorrect cuff size could result in erroneous and misleading BP measurement results

APPLICATION:
· Feed tab through the ring so that the artery label is visible on the outside of the cuff.

· Slide the cuff up the arm and align the artery marker over the brachial artery. Ensure that the ARTERY arrow is pointing down the arm.

· Cinch the cuff snugly on the arm and velcro closed. Tube should be pointing down the arm.

CLEANING PROCESS:

- Remove bladder.
- Soak in warm water for 5 minutes to loosen soil.
- Rinse under running warm water (50-140°F or 10-60°C) for 1 minute. A soft bristle brush may be used to remove visible soil.
- Fold up cuff so that hook is attached to loop with approximately 1/2" overlap and place in washing machine.
- Use a fresh solution of STERIS® Polystyca® 2X Concentrate Neutral Detergent.
- Machine wash the shell in warm water (50-140°F or 10-60°C).
- Remove shell from washing machine and rinse and massage under warm water for 5 minutes to remove any residual detergent. Run rinse water over both sides of the cuff as well as directly into the inner portion where the bladder is placed.
- Lay flat or line dry the shell.
- The bladder needs to be inserted back into the cuff sleeve so the pneumatic hose portion of the bladder is outside the sleeve.

NOTE: Do not place the cuff (sleeve) in a dryer.

DISINFECTION PROCESS:

NOTE: Before disinfecting, ensure that no liquid enters tubing by using a plug or taping off.

- Spray cuff with Quaternary Ammonium until soaked.
- Leave to soak for at least 10 minutes.
- Rinse with distilled water, ensuring liquid does not enter tube connector, and line dry.

Disposal: Please dispose of cuff according to local regulations

ENVIRONMENTAL CONDITIONS WARRANTY
SunTech Medical warrants our blood pressure cuff products to be free from defects in operating Ranges: 0 - 50°C Storage Ranges: -20 - 65°C

RH: 15 - 95% non-condensing

material and workmanship 24 months from the original date of purchase. This limited warranty covers the no charge replacement of the cuff under normal wear and tear conditions when returned to the attention of Service Department at the address below depending on location.

Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Service Department at the addresses below.

INTENDED USE & CLINICAL BENEFITS
The Eclipse, model SELF100, is intended to be used on pediatric and adult patients by trained clinicians, physicians, or a nurse. This cuff can be self-applied if user is trained by said professional. The cuff is not intended to be used on neonates. The clinical benefit of this cuff is its D-ring and pull-tab design specifically for self-application by the patient. This cuff is designed to stay in place during use to ensure consistent and accurate non-invasive blood pressure measurements. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the user's local competent authority.

CONTRAINDICATIONS
Contraindicated for use on neonates. Contraindicated for use with an MRI.

CE **MD**

SUNTECH
www.SunTechMed.com
SunTech Medical, Inc.
507 Airport Blvd., Suite 117
Morrissville, NC 27560, USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

EMERGO Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

UK **REP**
International Associates Limited
Centrum House, 38 Queen Street, Glasgow,
Lanarkshire, G1 3DX, UK
UKRP@ia-uk.com

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