
ENGLISH

Suntech Medical's All Purpose Clinical Grade Blood Pressure Cuffs

INDEX Line

Cuff index line must fall within range markings

Artery symbol and arrow should be placed over brachial or femoral artery

Symbol indicating product contains no PVC (polyvinyl chloride) material

Symbol indicating product contains no latex material

Product in compliance with Council Directive 93/42/EEC Medical Devices Directive

Symbol indicating arm circumference

Symbol indicating manufacturer

Symbol indicating lot code of cuff

Symbol indicating Caution: read the manual before use

RANGES / COLORS

Size:	Color:	Range:
Infant	Orange	8-13 cm
Child	Green	12-19 cm
Child Long	Green	12-19 cm
Sm Adult	Royal Blue	17-25 cm
Sm Adult Long	Royal Blue	17-25 cm
Adult	Navy Blue	23-33 cm
Adult Long	Navy Blue	23-33 cm
Lg Adult	Burgundy	31-40 cm
Lg Adult Long	Burgundy	31-40 cm
Thigh	Brown	38-50 cm

WARNINGS & INDICATIONS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb.

After washing ensure the size indication on the bladder and cuff shell match. Make sure that the cuff hose is threaded through one of the hose openings in the cuff.

Replace cuff when worn or damaged.

Do not machine wash bladder.

If patient has a mastectomy, place the cuff on the arm opposite the side of the mastectomy. If a patient has a double mastectomy, the cuff must be placed on the least dominate arm.

If a patient is undergoing intravascular access or therapy, or an arterio-venous (A-V) shunt, due to the presence of temporary interference of blood flow, do not apply the cuff on that arm. Potential injury to patient may occur.

The use of cuff while using other monitoring medical equipment could affect the accuracy or operation of said equipment.

Federal (US) law restricts this device to sale by or on the order of a physician.

Cuffs are intended for use with non-invasive blood pressure monitors only.

Follow the application instructions for use to ensure the correct size cuff for the patient. Failure to do so will adversely affect the accuracy of the reading.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Minimize limb movement during the measurement.

Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress. Do not compress or pinch cuff tubing.

Ensure patient's skin is intact before applying the cuff. Promptly remove cuff from patient when monitoring is not in progress.

During frequent or extended measurement, check cuff site and limb to ensure proper blood flow. Do not apply cuff where circulation may be compromised.

APPLICATION

- Place open cuff to the inner portion of the upper arm (or thigh) over the brachial artery.
- Align artery symbol ARTERY to the brachial (or femoral) artery.
- Use the RANGE indicator with the INDEX line to check that the arm falls within the specified range of that cuff. If it does not, select a cuff that better accommodates the limb circumference.
- Wrap the cuff snugly around the arm (or thigh).

Environmental Conditions

The following cleaning methods have been applied 20 times to the cuff without any apparent negative effects.

The cuff may be sprayed with a mild disinfectant solution (e.g. Cidezyme® ENZOL®, or 10% bleach solution), rinsed with distilled water and line dry. Ensure that no liquid enters tubing.

OR

To machine wash the cuff, remove the bladder and fully engage the hook and the loop. Machine wash warm with a mild detergent (50 – 130° F or 1 - 54° C) and line dry.

ENVIRONMENTAL CONDITIONS

Operating Ranges

Temp: 0 - 50°C

RH: 15 – 95% non-condensing

Storage Ranges

Temp: -20 - 65°C

RH: 15 – 95% non-condensing

WARRANTY

Suntech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. This limited warranty covers the no charge replacement of the cuff under normal wear and tear conditions when returned to the attention of Service Department at the address below depending on location.

Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Service Department at of the addresses below.

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SunTech Medical, Inc.
SunTech Medical, (Shenzhen) Co., Ltd.