# One Piece

## Durable BP Cuff

English

### Suntech Medical's One Piece Clinical Grade Blood Pressure Cuffs

- Index Line
- Cuff index line must fall within range markings
- Artery symbol and arrow should be placed over brachial or femoral artery
- Not made with PVC
- Not made with natural rubber latex
- Product in compliance with Council Directive 93/42/EEC Medical Devices Directive
- Symbol indicating arm circumference
- Symbol indicating lot code of cuff
- Symbol indicating manufacturer
- Symbol indicating Caution
- Indicates device sale only by or on the order of a licensed practitioner
- Indicates keep dry
- Medical Device

RANGES / COLORS		
Size:	Color:	Range:
Infant	Orange	8-13 cm
Child	Green	12-19 cm
Sm Adult	Royal Blue	17-25 cm
Sm Adult Plus	Royal Blue	18-29 cm
Adult	Navy Blue	23-33 cm
Adult Long	Navy Blue	23-33 cm
Adult Plus	Navy Blue	28-40 cm
Lg Adult	Burgundy	31-40 cm
Lg Adult Long	Burgundy	31-40 cm
Lg Adult Plus	Burgundy	40-55 cm
Thigh	Brown	38-50 cm

### WARNINGS

Pressurization of the cuff can temporarily cause loss of function if simultaneously used with monitoring equipment on the same limb. The cuff should not be placed on the patient's arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.

If luer lock connectors are used in the construction of tubing, there is a possibility to connect the cuff to intravenous fluids, allowing air to be pumped into a blood vessel, potentially causing serious injury.

Do not attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.

Following the application of the BP cuff, petechia formation (a minute reddish or purplish spot containing blood that appears in the skin) or Rumpel-Leede phenomenon (multiple petechia) on the arm, which may lead to idiopathic-thrombocytopenia (spontaneous persistent decrease in the number of platelets, associated with hemorrhagic conditions) or phlebitis (inflammation of a vein) may be observed.

As this medical device uses an alternative small-bore connector design different from those specified in the ISO 80369 series, there is a possibility that a misconnection can occur between this medical device and a medical device using a different alternative small-bore connector, which can result in a hazardous situation causing harm to the patient. Special measures need be taken by the user to mitigate these reasonable foreseeable risks.

### CAUTIONS

Federal (US) law restricts this device to sale by or on the order of a licensed healthcare practitioner. The cuff should not be applied over a wound as this can cause further injury.

Do not use cuff if there are any signs of damage. Failure to do so could affect measurement accuracy.

Minimize limb movement during the measurement.

Avoid contact with the cuff, other than that of the patient's limb, while measurement is in progress.

A compressed or kinked connection hose may cause continuous cuff pressure resulting in blood flow interference and potentially harmful injury to the patient.

Using an incorrect cuff size could result in erroneous and misleading BP measurement results.

### APPLICATION

Follow the application instructions for use to ensure the correct size cuff for the patient. Failure to do so will adversely affect the accuracy of the reading.

- Place open cuff around the inner portion of the upper arm (or thigh).
- Align artery symbol **ARTERY** to the brachial (or femoral) artery.
- Use the **RANGE** indicator with the **INDEX** line to check that the limb falls within the specified range of that cuff. If it does not, select a cuff that better accommodates the limb circumference.
- Wrap the cuff snugly around the arm (or thigh).

Promptly remove cuff from patient when monitoring is not in progress.

### CLEANING INSTRUCTIONS

The following cleaning methods have been applied 20 times to the cuff without any apparent negative effects.

The cuff may be sprayed with a mild disinfectant solution (e.g. Cidezyme®, ENZOL®, or 10% bleach solution), rinsed with distilled water and line dry. Ensure that no liquid enters tubing.

OR  
The cuff may be wiped down with a mild disinfectant wipe (e.g. Sani-Cloth®) and line dry.

Disposal: Please dispose of cuff according to local regulations.

### ENVIRONMENTAL CONDITIONS

	Operating Ranges
	Temp: 0 - 50°C RH: 15 – 95% non-condensing
	Storage Ranges
	Temp: -20 - 65°C RH: 15 – 95% non-condensing

### WARRANTY

Suntech Medical warrants our blood pressure cuff products to be free from defects in material and workmanship 24 months from the original date of purchase. This limited warranty covers the no charge replacement of the cuff under normal wear and tear conditions when returned to the attention of Service Department at the address below depend-ing on location. Contact the Service Department at either of these locations to receive a Return Material Authorization number before sending any product. Cuffs should be returned to the attention of the appropriate Ser-vice Department at of the addresses below.

### OPD CUFF INTENDED USE STATEMENT

The One Piece Durable Cuff is intended to be used with a manual or automatic non-invasive sphygmomanometer to determine blood pressure parameters on pediatric and adult patients. They are not intended to be used on neonates. They are intended to be used with prescription (clinical grade) blood pressure monitors; they also could be used with Over-The-Counter (home) blood pressure monitors.

Suntech Medical, Inc.  
507 Airport Boulevard, Suite 117  
Morrisville, NC 27560-8200 USA  
Tel: + 1.919.654.2300  
1.800.421.8626  
Fax: + 1.919.654.2301

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