

Soft, Single Patient-Use
Blood Pressure Cuffs


SunTech® *Soft Disposable*
SINGLE-PATIENT USE BP CUFF

Providing Affordable Infection Control



Tapered End & Rounded Edges



External & Internal Range Indicators



Patient ID for Single Patient-Use



- + Engineered for patient comfort
- + 100% leak tested

- + Cuff Standardization Program
- + Reliability from the BP experts

For over twenty years, SunTech Medical® has focused solely on the development of blood pressure measurement technologies for the most challenging and complex clinical environments. As automated monitors have become routine for blood pressure measurement, it is important to remember that the cuff is the sensor. With this in mind, we have applied our expertise to engineer a new line of disposable blood pressure cuffs. These cuffs are designed with a tapered end to improve patient comfort and external and internal range indicators to ensure accuracy and proper fit. This unique approach makes SunTech Medical disposable cuffs the easy choice to meet your infection control requirements at an exceptional value.



Full Range of Cuff Sizes, Including Neonatal, Infant & Child



One or Two-Tube Designs Available








+ Features & Benefits





- Patient Comfort:**
- Economical:**
- Ease of Use:**
- Patient Identification:**
- Reliable Performance:**
- Regulatory Requirements:**
- Durability:**
- Latex-Free:**
- Cuff Standardization Program (CSP):**

+ Description







- Softer materials, tapered end and rounded edges provide a gentle cuff solution.
- Priced to support single-patient use, particularly where infection control is a concern.
- Tapered end for easier cuff application. Color-coded for quick size selection. Interior & exterior range indicators ensure accurate cuff sizing and fit.
- Dedicated location for single patient identification. Minimizes cross-contamination.
- 100% leak tested for trust, value, and reliability.
- Engineered to meet AHA, AAMI-SP10, and MDD guidelines for results you can trust.
- Value in a durable, single-piece design for multiple inflations.
- 100% latex-free materials promotes patient safety.
- Expertise in standardizing your facility to one BP cuff connector: one cuff, any patient, any monitor, any department.

+ Ranges

-  Neonatal 1
3-6 cm
-  Neonatal 2
4-8 cm
-  Neonatal 3
6-11 cm
-  Neonatal 4
7-13 cm
-  Neonatal 5
8-15 cm
-  Infant
8-13 cm
-  Child/Child Long
12-19 cm

-  Small Adult/
Small Adult Long
17-25 cm
-  Adult/Adult Long
23-33 cm
-  Large Adult/
Large Adult Long
31-40 cm
-  Thigh
38-50 cm

+ Symbols

-  Index Line.
-  Cuff index line must fall within range markings for an accurate measurement.
-  Artery symbol and arrow should be placed over brachial or femoral artery.
-  Product contains no latex material.
-  Product in compliance with Council Directive 93/42/EEC Medical Devices Directive.
-  Arm circumference



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235

www.SunTechMed.com

SunTech Medical, Ltd.
2/F of Building A, Jinxiongda Technology Park
Guanlan, Bao'an District Shenzhen, 518110, PRC
Tel: +86-755-2958 8810
+86-755-29588986 (Sales)
+86-755-29588665 (Service)
Fax: +86-755-2958 8829

www.CN.SunTechMed.com

ETL CLASSIFIED



Intertek
8785
MEDICAL
ELECTRICAL
EQUIPMENT



0413