


Technical Data Sheet

All-Purpose Cuff DURABLE TWO-PIECE



Cuff Size	1-Tube or 2-Tube	Range (cm)	Bladder Length x Width (mm)	Overall Length x Width (mm)	Female Quick Connect	Male Quick Connect	Threaded (Screw)	Male Bayonet (Metal)	Male Bayonet (Plastic)	Female & Male Quick Connects	No Connector
Infant	1-Tube	8-13	110 x 38	230 x 63	98-0080-01	98-0084-31	98-0084-21	98-0084-41	98-0084-71	N/A	98-0084-11
Infant	2-Tube	8-13	110 x 38	230 x 63	98-0082-01	98-0085-31	98-0085-21	98-0085-41	98-0085-71	98-0058-61	98-0085-11
Child	1-Tube	12-19	150 x 68	315 x 93	98-0080-02	98-0084-32	98-0084-22	98-0084-42	98-0084-72	N/A	98-0084-12
Child	2-Tube	12-19	150 x 68	315 x 93	98-0082-02	98-0085-32	98-0085-22	98-0085-42	98-0085-72	98-0058-62	98-0085-12
Child Long	1-Tube	12-19	150 x 68	445 x 93	98-0080-03	98-0084-33	98-0084-23	98-0084-43	98-0084-73	N/A	98-0084-13
Child Long	2-Tube	12-19	150 x 68	445 x 93	98-0082-03	98-0085-33	98-0085-23	98-0085-43	98-0085-73	98-0085-63	98-0085-13
Small Adult	1-Tube	17-25	193 x 93	440 x 118	98-0080-04	98-0084-34	98-0084-24	98-0084-44	98-0084-74	N/A	98-0084-14
Small Adult	2-Tube	17-25	193 x 93	440 x 118	98-0082-04	98-0085-34	98-0085-24	98-0085-44	98-0085-74	98-0085-64	98-0085-14
Small Adult Long	1-Tube	17-25	193 x 93	515 x 118	98-0080-05	98-0084-35	98-0084-25	98-0084-45	98-0084-75	N/A	98-0084-15
Small Adult Long	2-Tube	17-25	193 x 93	515 x 118	98-0082-06	98-0085-35	98-0085-25	98-0085-45	98-0085-75	98-0085-65	98-0085-15
Adult	1-Tube	23-33	265 x 125	505 x 150	98-0080-06	98-0084-36	98-0084-26	98-0084-46	98-0084-76	N/A	98-0084-16
Adult	2-Tube	23-33	265 x 125	505 x 150	98-0082-07	98-0085-36	98-0085-26	98-0085-46	98-0085-76	98-0085-66	98-0085-16
Adult Long	1-Tube	23-33	265 x 125	645 x 150	98-0080-07	98-0084-37	98-0084-27	98-0084-47	98-0084-77	N/A	98-0084-17
Adult Long	2-Tube	23-33	265 x 125	645 x 150	98-0082-08	98-0085-37	98-0085-27	98-0085-47	98-0085-77	98-0085-67	98-0085-17
Large Adult	1-Tube	31-40	320 x 158	645 x 183	98-0080-08	98-0084-38	98-0084-28	98-0084-48	98-0084-78	N/A	98-0084-18
Large Adult	2-Tube	31-40	320 x 158	645 x 183	98-0082-09	98-0085-38	98-0085-28	98-0085-48	98-0085-78	98-0085-68	98-0085-18
Large Adult Long	1-Tube	31-40	320 x 158	810 x 183	98-0080-09	98-0084-39	98-0084-29	98-0084-49	98-0084-79	N/A	98-0084-19
Large Adult Long	2-Tube	31-40	320 x 158	810 x 183	98-0082-09	98-0085-39	98-0085-29	98-0085-49	98-0085-79	98-0085-69	98-0085-19
Thigh	1-Tube	38-50	390 x 195	760 x 220	98-0080-10	98-0084-40	98-0084-30	98-0084-50	98-0084-80	N/A	98-0084-20
Thigh	2-Tube	38-50	390 x 195	760 x 220	98-0082-10	98-0085-40	98-0085-30	98-0085-50	98-0085-80	98-0085-70	98-0085-20

One-Piece Cuff

DURABLE BLADDERLESS



Cuff Size	1-Tube or 2-Tube	Range (cm)	Bladder Length x Width (mm)	Overall Length x Width (mm)	Female Quick Connect	Threaded	Bayonet	No Connector
Infant	1-Tube	8-13	120 x 54	218 x 54	98-0600-00	98-0600-20	98-0600-40	98-0600-50
Infant	2-Tube	8-13	120 x 54	218 x 54	98-0600-10	98-0600-30	N/A	98-0600-60
Child	1-Tube	12-19	160 x 83	304 x 83	98-0600-01	98-0600-21	98-0600-41	98-0600-51
Child	2-Tube	12-19	160 x 83	304 x 83	98-0600-11	98-0600-31	N/A	98-0600-61
Small Adult	1-Tube	17-25	203 x 108	430 x 108	98-0600-03	98-0600-23	98-0600-43	98-0600-53
Small Adult	2-Tube	17-25	203 x 108	430 x 108	98-0600-13	98-0600-33	N/A	98-0600-63
Adult	1-Tube	23-33	275 x 140	495 x 140	98-0600-05	98-0600-25	98-0600-45	98-0600-55
Adult	2-Tube	23-33	275 x 140	495 x 140	98-0600-15	98-0600-35	N/A	98-0600-65
Adult Long	1-Tube	23-33	275 x 140	630 x 140	98-0600-06	98-0600-26	98-0600-46	98-0600-56
Adult Long	2-Tube	23-33	275 x 140	630 x 140	98-0600-16	98-0600-36	N/A	98-0600-66
Large Adult	1-Tube	31-40	330 x 173	635 x 173	98-0600-07	98-0600-27	98-0600-47	98-0600-57
Large Adult	2-Tube	31-40	330 x 173	635 x 173	98-0600-17	98-0600-37	N/A	98-0600-67
Large Adult Long	1-Tube	31-40	330 x 173	805 x 173	98-0600-08	98-0600-28	98-0600-48	98-0600-58
Large Adult Long	2-Tube	31-40	330 x 173	805 x 173	98-0600-18	98-0600-38	N/A	98-0600-68
Thigh	1-Tube	38-50	400 x 210	745 x 210	98-0600-09	98-0600-29	98-0600-49	98-0600-59
Thigh	2-Tube	38-50	400 x 210	745 x 210	98-0600-19	98-0600-39	N/A	98-0600-69

Cuff Component	Material	All-Purpose	One-Piece
Cuff Shell	Nylon, 200 Denier	•	•
Hook	Nylon	•	•
Loop	Nylon	•	•
Thread	Nylon	•	
Bladder	Nylon, 200 Denier	•	
Bladder Tubing	TPE (thermoplastic elastomer)	•	•
Bladder Nipple	Polyurethane	•	•
Outer Label	Polyester	•	
100% Latex-free		•	•
100% PVC-free		•	•
100% DEHP-free		•	•
Regulatory Requirements			
Meets AHA, AAMI-SP10 Guidelines		•	•
In compliance with Council Directive 93/42/EEC Medical Devices Directive		•	•
Meets ISO 10993 standard		•	•



<p>SunTech Medical, Inc. 507 Airport Boulevard, Suite 117 Morrisville, NC 27560-8200 USA Tel: + 1.919.654.2300 1.800.421.8626 Fax: + 1.919.654.2301</p> <p>www.SunTechMed.com</p>	<p>SunTech Medical, Ltd. Oakfield Estate Eynsham, Oxfordshire England Tel: + 44 (0) 1865.884.234 Fax: + 44 (0) 1865.884.235</p> <p>www.SunTechMed.com</p>	<p>SunTech Medical, Ltd. 2/F of Building A, Jinxiongda Technology Park Guanlan, Bao'an District Shenzhen, 518110, PRC Tel: +86-755-2958 8810 +86-755-29588986 (Sales) +86-755-29588665 (Service) Fax: +86-755-2958 8829</p> <p>www.SunTechMed.com/CN</p>
---	---	--