



Introducing

Oscar 2TM

**SphygmoCor
inside**



Oscar 2
SphygmoCor
inside



24-Hour Ambulatory Blood Pressure Monitoring System
with Pulse Wave Analysis

- + Backed by 800 Published Studies
- + Fully-Featured Physiological Waveforms

- + AccuWin ProTM 4 – ABPM Analysis & Reporting
- + Industry-Leading Measurement Performance

+ SphygmoCor Inside

The gold standard in Pulse Wave Analysis & Central Blood Pressure Measurement

- Validated to inter-arterial catheter
- Over 800 published studies
- Sub-Diastolic CBP Waveform Capture
- Fully featured central blood pressure waveforms



+ Advanced Analysis with AccuWin Pro™ 4

- Automated interpretation of ABPM data
- Programmable quality control metrics
- Detailed pulse wave analysis (PWA)
- HL7 Output Capability
- Automated White-coat analysis

+ Specifications

Validations	BHS (A/A), ESH, ANSI/AAMI 81060-2
BP Technique	Oscillometry with step deflation
BP Range	25-260 mmHg
HR Range	40-200 bpm
Sampling Periods	3 adjustable (awake, asleep, and special)
Time Intervals	5, 10, 15, 20, 30, 45, 60, 90 and 120 min.
Event Key	Program to allow additional patient readings
Data Storage	Flash memory stores up to 250 readings
Power Source	Two AA batteries
Weight	275 grams, including batteries
Dimensions	120 x 70 x 30 mm
PC Requirements	USB port, Windows® XP or later

+ SunTech® ABPM Blood Pressure Cuffs

- Patented stretch sleeve keeps the cuff in place both day and night
- Optimal patient comfort and data validity
- 4 cuff sizes to fit almost any arm (18cm-55cm)
- Removable bladder for machine washing

+ The Oscar 2™ with SphygmoCor Inside System Includes

- Oscar 2™ ABP monitor
- AccuWin Pro™ 4 software
- 2 SunTech® ABPM cuffs (Sizes 2 and 3)
- Carrying pouch & patient belt
- USB data cable
- 4 AA alkaline batteries
- User's Manual & Quick Start Guide

Part Numbers	Descriptions
99-0133-03	Oscar 2™ with SphygmoCor Inside ABPM System
98-0239-01	ABPM Cuff Size 1 (18-26 cm)
98-0239-02	ABPM Cuff Size 2 (25-35 cm)
98-0239-03	ABPM Cuff Size 3 (32-44 cm)
98-0239-04	ABPM Cuff Size 4 (42-55 cm)



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235



www.SunTechMed.com

SunTech Medical, Ltd.
105 HuanGuan South Road, Suite 15 2~3/F
DaHe Community Guanlan,
LongHua District, Shenzhen
GuangDong PRC 518110
Tel: + 86-755-2958 8810
+ 86-755-29588986 (Sales)
+ 86-755-29588665 (Service)
Fax: + 86-755-2958 8829
www.SunTechMed.com.CN

ETL CLASSIFIED



Intertek
1795
MEDICAL
ELECTRICAL
EQUIPMENT



0413

82-0404-00-MA-RevA