

IHE Integration Statement			
Vendor	Product Name	Version	Date
SunTech Medical, Inc.	MODEL 260	1.0	3/16/16

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Consistent Time	Time Client	None
Device Enterprise Communication	Device Observation Reporter	None
Patient Demographics Query	Patient Demographics Consumer	None

Internet address for vendors IHE information : http://www.suntechmed.com	
Links to Standards Conformance Statements for the Implementation	
HL7	
DICOM	
Link to general information on IHE : www.ihe.net	