



 SunTech® **Vet20™**
BP MONITOR



Spot-Check Blood Pressure Monitor
for Companion Animals

SunTech Vet20

The portable SunTech® Vet20™ quickly and easily provides reliable blood pressure (BP) measurements on awake companion animals, ensuring a simpler assessment process for veterinarians and their staff members. Utilizing SunTech's trusted Advantage™ VET BP technology, the Vet20 quietly performs motion-tolerant, oscillometric BP measurements, delivering accurate results. The intuitive touch screen interface includes all of the necessary features for reviewing and analyzing BP data. Alarms can be silenced to avoid frightening the animal and saved measurement data is always retained when the device is powered off.

Standard Features

- Reliable:** As a result of SunTech's motion-tolerant technology, the Vet20 has a high success rate for BP readings.
- Simple to Use:** Intuitive, touch screen interface has all necessary features for reviewing and analyzing BP data.
- Accurate:** BP algorithm was specifically developed for companion animals by experts at a leading vet school.
- Silent:** Device quietly measures BP and alarms may be silenced to avoid frightening the animal.
- Fast:** Take multiple measurements in a matter of minutes!
- Portable:** Battery-powered device can be easily moved between exam rooms.
- Memory:** The Vet20's advanced memory ensures data is always retained when the device is powered off.
- Averaging:** Users can select specific measurements from the device memory and perform automated BP averaging.



Technical Specifications

The Vet20 Includes:

- Vet20 System Monitor
- Protective Armour in Flamingo Pink, Peacock Blue, or Tree Frog Green
- Veterinary Blood Pressure Cuffs Sizes 1-6
- Patient Hose, 4 ft (1.2 m)
- AC Adapter

General Specifications

Systolic:	40 to 265 mmHg
Diastolic:	20 to 200 mmHg
Mean Arterial Pressure:	27 to 222 mmHg
Heart Rate:	25 to 300 BPM
Device Weight:	1.3 lbs (with battery), 590 grams (with battery)
Dimensions:	6.25" x 5" x 5.25", 15.9 cm x 12.7 cm x 13.3 cm
Warranty:	2 years on monitor 90 days on cuffs/hose

Accessories

SunTech Vet20 Armour

The removable protective armour is available in 3 colors and features a convenient handle.



Flamingo Pink



Peacock Blue



Tree Frog Green

AC Adapter

While the portability of the battery-powered Vet20 is convenient, the AC Adapter accessory is a great power option when your Vet20 stays in one place. The A/C Adapter may be used with any SunTech Vet20 BP monitor and includes 4 plug configurations to meet most international power requirements.

SunTech Veterinary Cuffs



Soft materials and rounded corners provide a gentle cuff solution. SunTech Medical BP cuffs are color-coded for quick size selection and range indicators ensure accurate cuff sizing and fit. Included with the Vet20 BP device are cuffs size 1–6. **SunTech Vet BP cuffs can be used with other veterinary BP monitors! Ask about our SunTech Vet Cuff Conversion Kit.**


SunTech Medical®
The Difference in Clinical Grade®

With more than 30 years dedicated to developing our world-class blood pressure technology, we are confident that we don't miss a beat - so you won't either. It's our legacy. It's SunTech Medical. ***It's The Difference In Clinical Grade.***



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: +1.919.654.2301
SunTechMed.com

EC REP

SunTech Medical, Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235
SunTechMed.com

SunTech Medical (Shenzhen) Co., Ltd.
105 HuanGuan South Road, Suite 15 2~3/F
DaHe Community Guanlan,
LongHua District, Shenzhen
GuangDong PRC 518110
Tel.: + 86.755.29588810
+ 86.755.29588986 (Sales)
+ 86.755.29588665 (Service)
Fax: + 86.755.29588829
SunTechMed.com.CN





SunTech Vet20 Part Numbers

Vet20 System Monitors

99-0135-06	Vet20 Monitor, Flamingo Pink Armour
99-0135-07	Vet20 Monitor, Peacock Blue Armour
99-0135-12	Vet20 Monitor, Tree Frog Green Armour

SunTech Veterinary BP Cuffs

98-0400-80	Cuff #1	Box of 20	Range 3-6 cm
98-0400-81	Cuff #2	Box of 20	Range 4-8 cm
98-0400-82	Cuff #3	Box of 20	Range 6-11 cm
98-0400-83	Cuff #4	Box of 20	Range 7-13 cm
98-0400-84	Cuff #5	Box of 20	Range 8-15 cm
98-0400-F1	Cuff #6	Box of 20	Range 12-19 cm
98-0400-F3	Cuff #7	Box of 20	Range 17-25 cm
98-0240-00	Cuffs 1-6	Pack of 6	
91-0028-72	Patient Hose, 4 ft (1.2 m)		
91-0028-75	Patient Hose, 6 ft (1.8 m)		
91-0028-76	Patient Hose, 10 ft (3 m)		

Accessories

39-0195-03	Flamingo Pink Armour
39-0195-04	Peacock Blue Armour
39-0195-02	Tree Frog Green Armour
19-0018-00	AC Adapter

Extended Warranty Program

83-0021-00	One Year Extended Warranty
83-0021-01	Three Years Extended Warranty

Extended Warranty Terms and Conditions

The Extended Warranty can be purchased either:

- At the time of the original sale
- Prior to the end of the initial warranty period