

Memo

To: SunTech Customers
From: SunTech Medical Sales Team
Date: 9/8/2021
Re: Oscar 2 Classic and Orbit Cuffs End of Life - IMPORTANT UPDATE

Dear SunTech Customer,

The SunTech Oscar 2 Classic (Model 222) went End of Life on December 16, 2018. The Oscar 2 (Model 250) is available as a replacement product and is backwards compatible with AccuWinPro v3. Affected Oscar 2 Classic Part Numbers are any part numbers in the range 99-0012-00 to 99-0012-99.

The SunTech line of Classic Orbit cuffs is no longer available for purchase with immediate effect. The newest Orbit cuff range is compatible with all Oscar 2 models and can be used as replacement cuffs. The new Orbit cuff range uses the latest AHA guidelines on cuff size and uses a numeral size system. Affected Orbit Cuff Part Numbers are any part numbers in the ranges 98-0063-00 to 98-0063-99 and any part numbers in the ranges 98-0122-00 to 98-0125-99

The following replacement Orbit cuff part numbers should be used in place of Orbit Classic cuffs:

Classic Orbit Cuff			New Replacement Orbit Cuff		
Size	Range	Part Number	Size	Range	Part Number
Small Adult	18cm – 27cm	98-0063-11	Size 1	18cm – 26cm	98-0239-01
Adult	25cm – 35cm	98-0063-12	Size 2	25cm – 35cm	98-0239-02
Adult Plus	33cm – 40cm	98-0063-13	Size 3	32cm – 44cm	98-0239-03
Large Adult	39cm – 46cm	98-0063-14	Size 4	42cm – 55cm	98-0239-04

More information on the newest Orbit cuff is available on our [website](#).

UPDATE: Repair, support and replacement parts previously were stated to remain available for this product for a period of 7 years from sale. Due to market conditions beyond our control, we can no longer guarantee availability of replacement parts.

For questions on availability and end of life, contact Rod Thomson at rthomson@suntechmed.com

Regards,
SunTech Team

US Phone Contact: +1.919.654.2300 (US)
+1.800.421.8626 (US Toll Free)

Email Contact: rthomson@SunTechMed.com