

Short-Term Use
Blood Pressure Cuffs

SunTech® Vinyl

SHORT-TERM USE BP CUFFS

Durable and Economical Enough for Short-Term Use



Full Range of Cuff Sizes



External & Internal Range Indicators



Tapered End & Rounded Edges



- + Easily Cleaned & Disinfected
- + 100% latex-free

- + 100% leak tested
- + Engineered for results you can trust

For nearly thirty years, SunTech Medical® has focused solely on the development of blood pressure measurement technologies for the most challenging and complex clinical environments. As automated monitors have become routine for blood pressure measurement, it is important to remember that the cuff is the sensor. With this in mind, SunTech now offers an economical blood pressure cuff solution for clinical areas where cleaning and disinfection are important—during procedures, in the emergency department, and EMS. SunTech's Vinyl BP Cuffs can be easily cleaned, are rugged enough for short-term use, and are less expensive than traditional durable cuffs.

+ Features & Benefits

Cleanable:

Economical:

Durability:

Ease of Use:

Reliable Performance:

Regulatory:

Latex-Free:

+ Description

Cuffs are made with vinyl material that can be easily cleaned and disinfected.

Lower cost than traditional durable cuffs makes short-term use affordable.

Value in a durable, single-piece design for multiple inflations.
















Tapered end for easier cuff application. Color-coded for quick size selection. Interior & exterior range indicators ensure accurate cuff sizing and fit.

100% leak tested for trust, value, and reliability.

Engineered to meet AHA, AAMI-SP10, and MDD guidelines for results you can trust.

100% latex-free materials promote patient safety.

+ Ranges & Design Specifications

+ Ranges & Design Specifications					One Tube							Two Tube	
													
Range Symbols	Cuff Size (Print Color)	Range (cm)	Bladder Length x Width (mm)	Overall Length x Width (mm)	Female Quick Connect	Threaded (Screw)	Male Bayonet (Plastic)	No Connector	Female Quick Connect	Male Bayonet (Plastic)	No Connector		
	Infant (Orange)	8-13	120 x 54	228 x 59	98-0800-10	98-0800-00	98-0800-30	98-0800-50	98-0800-20	98-0800-40	98-0800-60		
	Child (Green)	12-19	160 x 83	314 x 88	98-0800-11	98-0800-01	98-0800-31	98-0800-51	98-0800-21	98-0800-41	98-0800-61		
	Small Adult (Turquoise)	17-25	203 x 108	440 x 113	98-0800-13	98-0800-03	98-0800-33	98-0800-53	98-0800-23	98-0800-43	98-0800-63		
	Adult (Navy Blue)	23-33	275 x 140	505 x 145	98-0800-15	98-0800-05	98-0800-35	98-0800-55	98-0800-25	98-0800-45	98-0800-65		
	Adult Long (Navy Blue)	23-33	275 x 140	640 x 145	98-0800-16	98-0800-06	98-0800-36	98-0800-56	98-0800-26	98-0800-46	98-0800-66		
	Large Adult (Burgundy)	31-40	330 x 173	645 x 178	98-0800-17	98-0800-07	98-0800-37	98-0800-57	98-0800-27	98-0800-47	98-0800-67		
	Large Adult Long (Burgundy)	31-40	330 x 173	815 x 178	98-0800-18	98-0800-08	98-0800-38	98-0800-58	98-0800-28	98-0800-48	98-0800-68		
	Thigh (Brown)	38-50	400 x 210	755 x 215	98-0800-19	98-0800-09	98-0800-39	98-0800-59	98-0800-29	98-0800-49	98-0800-69		



SunTech Medical, Inc.
507 Airport Boulevard, Suite 117
Morrisville, NC 27560-8200 USA
Tel: + 1.919.654.2300
1.800.421.8626
Fax: + 1.919.654.2301

www.SunTechMed.com

SunTech Medical, Ltd.
Oakfield Industrial Estate
Stanton Harcourt Road
Eynsham, Oxfordshire OX29 4TS
England
Tel: + 44 (0) 1865.884.234
Fax: + 44 (0) 1865.884.235
www.SunTechMed.com

SunTech Medical, Ltd.
2/F of Building A, Jinxiongda Technology Park
Guanlan, Bao'an District Shenzhen, 518110, PRC
Tel: + 86-755-2958 8810
+ 86-755-29588986 (Sales)
+ 86-755-29588665 (Service)
Fax: + 86-755-2958 8829
www.SunTechMed.com.cn

ETL CLASSIFIED



Intertek
7195
MEDICAL
ELECTRICAL
EQUIPMENT



0413
82-0130-01 Rev. A